

## Setting up Your System

Press and hold the **S** button for 3 seconds to enter the SET mode. After all settings are finished, press and hold the **S** button for 3 seconds to turn off the meter.

Press **>>** to reach the correct value. Press and hold **>>** to scroll faster.

### Adjusting the Date and Time

#### Step 1 Entering the SET Mode


Press and hold the **S** button for 3 seconds to enter the SET mode. After all the segments flash across the screen, the 'SET' character will be displayed on the screen.

Press the **S** button again to progress to setting the year.



#### Step 2 Setting the Year

Press and release **>>** to adjust until the correct year appears. After setting the year, press the **S** button to confirm your selection and progress to setting the month.



#### Step 3 Setting the Month

A number indicating the month will blink on the left corner of the screen. Press **>>** until the correct month appears. Press the **S** button to confirm your selection and progress to setting the date.



#### Step 4 Setting the Date

Press **>>** until the screen displays the correct date. Press the **S** button to confirm the date and progress to setting the time.



#### Step 5 Setting the Time Format

The meter can be set in the AM/PM 12-hour or the 24-hour clock format. Press **>>** to select a format. The AM•PM symbol is not displayed in the 24-hour format. After selecting the format press the **S** button to progress to setting the hour.



### Step 6 Setting the Hour

Press **>>** until the correct hour appears.  
After the hour is set, press the **S** button to progress to setting the minute.



### Step 7 Setting the Minute

Press **>>** until the correct minute appears.  
After setting the minute, press the **S** button to progress to setting the sound.



### Setting the Sound On/OFF

#### Step 8

On pressing **>>**, the screen will display On or OFF. Press **S** button to confirm the selection.




The meter will beep in the following instances, if set to On.

- When you push the **S** button or **>>** button to turn on the meter
- When the test strip is inserted in the meter
- When the blood sample is absorbed into the test strip and the test starts
- When the test result is displayed
- When you push the **>>** button to set the post-meal (PP2) alarm
- When it is time for a preset blood glucose test



If the sound is set to OFF, none of the sound functions will work.

After setting the sound, press the **S** button to progress to setting the strip expiration date indicator.


**Note:**  symbol appears only when the sound is set to OFF.

## Turning on the Strip Expiration Date Indicator

### Step 9

This mode allows you to turn the strip expiration date indicator on or off. This mode turns the function on or off only. See page 40 to set the strip expiration date.

When "EP" blinks on the screen, press **>>**. The screen will display "On" or "OFF". Press the **S** button to confirm the setting and to progress on to setting the hypoglycaemia indicator.




**Note:** If the pre-set expiration date expires, the meter will display the following. For example, in the case when the expiration date is set to October of 2014, the meter displays "EP" at the start of November of 2014.



## Turning on the Hypoglycemia (HYPO) Indicator


### Step 10

In this mode the Hypoglycemia (possible low blood sugar) level can be selected.

The 'HYPO' character will be displayed with hypoglycemia symbol . On pressing the **>>** button, the screen will display the On or OFF. Press the **S** button when 'On' appears to enter the setting.

Press the **>>** button until the correct hypoglycemia level between 1.1 and 5.0 mmol/L is set. Then, press the **S** button to confirm the level. If you do not want to set the indicator, press the **S** button while the screen displays 'no'. Then, the screen will return to step 2. See page 13.

To confirm all settings, press and hold the **S** button for 3 seconds and the screen will switch off.



**Note:** If the test result is lower than the pre-set hypoglycemia level, the meter will display the following.



**Caution:** Before the HYPo indicator is set, ask your doctor or diabetes nurse to help you decide what your hypoglycemia level will be.