

Te Arataki'anga Viviki

CareSens™ NPOP

Te Akanoonoo'anga no te
Akara Vaito i te Tuka i Roto i te Toto



PGJ5E3033 REV1 03/2014


4. Te teatea-mamao'anga i te Matīni Ko'i Toto Akara'anga

- Takaviri taratarā'ia ma te akapae i te ope te ka rauka i te akatuke'ia.
- Tāru piri ketaketa pu'ia tetai mātipi meangiti ou ki roto i te mouranga o te mātipi meangiti. Mou piri meitaki ia te mātipi meangiti. Takaviri marie ia no te 'u'uti akapae mai i te apinga paruru tapoki. Akono meitaki ia te apinga paruru tapoki no te ta-angaanga akaou anga i te ope o te mātipi meangiti. Akaoki'ia te ope o te matīni ko'i toto.
- Iki'ia te 'o'onu tei anoano'ia mei te ta'i ki te itu (1 – 7) o te ope te ka rauka i te akatuke'ia, o te matīni ko'i toto akara'anga. Iki'ia te 'o'onu na roto i te takavirivir'i'anga ia a runga i te ope te ka rauka i te akatuke'ia, kia āiteite te numero tei akano'ia ki te apinga to'u (arrow). Kua tamanako'ia kia akano'ia te akamata'anga ki te toru (3).
- Akapano'ia te matīni ko'i toto akara'anga na roto i te 'u'uti'anga i te paero pātekateka kia rongo'ia atu te tangitangi me tangitangi mai. Ka muramura te māramarama akatinamou'anga tei akaki'ia.

5. Te Teatea-mamao'anga i te Mita e te Potonga Apinga Akara Akatinamou'anga

- Pu'a tāmā'ia toou nga rima e te ngai o te apinga akara'anga. Tāmā meitaki ia ki te vai, ka tāmā meitaki tikai atu ei.
- Tāru'ia tetai potonga apinga akara akatinamou'anga ma te uri ki runga, ki roto i te au tu'anga arairai. Ka ata mai te akairo ▲ ki runga i te ngai akaariari.


1. Te Mita CareSens N POP no te Tuka i Roto i te Toto




2. Te Potonga Apinga Akara Akatinamou CareSens N POP




3. Te Matīni Ko'i Toto Akara'anga



6. Te Raverave Kapiti Atu Anga i te Toto



7. Te Akapae Titiri'anga i te Au Mātipi Meangiti tei Oti i te Ta-angaanga'ia



8. Te Akaari'anga e te Tuatua Akakitekite

- Er 1**: Kua tāru'ia tetai potonga o te apinga akara akatinamou'anga tei oti i te ta-angaanga'ia.
-> Tāoki akaou'ia te akara akatinamou'anga ki tetai potonga ou.
- Er 2**: Kua kapiti'ia atu te toto, me kore ra, te vairakau akateretere akara'anga, i mua ake i ata mai ei te akairo ▲.
-> Tāoki akaou'ia te akara akatinamou'anga ki tetai potonga ou ma te tiaki kia atu mai te akairo ▲ i mua ake ka kapiti atu ei i te toto me kore ra, te vairakau akateretere akara'anga.
- Er 4**: Kua teitei kino te turanga visokositi (viscosity) o te toto, me kore ra, kare e rava meitaki ra te turanga toto.
-> Tāoki akaou'ia te akara akatinamou'anga i muri ake ka tāru'ia ai tetai potonga ou.
- Er 5**: Penei, ka ata mai teia tuatua akakite me tarevake te potonga apinga akara akatinamou o te turanga tuka i roto i te toto, e kare te reira i te potonga apinga akara akatinamou toto a te CareSens N.
-> Tāoki akaou'ia te akara akatinamou'anga ki tetai potonga akara akatinamou CareSens N.
- Er 6**: E manamanan tei tupu ki te mita.
-> Auraka e ta-angaanga i te mita. Kapiki'ia atu to kotou tangata o te ngai angaanga oko o te i-SENS.
- HIC**: Tei runga atu te vāito reva i te tuatau i raveia ai te akara akatinamou'anga.
-> E neke ki tetai ngai tei tau te vāito reva ki tei matau'ia i te rave i teia (5 – 50°C), ma te tāoki akaou i te akara akatinamou'anga i muri ake me tae te vāito reva o te mita e te au potopotonga apinga, ki te turanga vāito reva tei tau no te rave i teia.
- LoC**: Tei raro mai te vāito reva i te tuatau i raveia ai te akara akatinamou'anga i te turanga tei matau'ia me raveia teia.

Te Arataki'anga Viviki

CareSens™ NPOP
Te Akanoonoo'anga no te
Akara Vāito i te Tuka i Roto i te Toto

